



Technical Publications

Direction 5215475-100

Revision 3

Reporting Tool 2.7

DICOM CONFORMANCE STATEMENT

Copyright © 2005-2010 by General Electric Co.

Do not duplicate



REVISION HISTORY

Revision	Date	Reason for Change
A	Sep 25, 2007	Initial revision
1	Oct 19, 2007	After review
2	Mar 06, 2009	<p>Added support for</p> <ul style="list-style-type: none"> • Mammography CAD SR IOD • Key Object Selection Document IOD • DICOM Encapsulated PDF IOD <p>Added changes after review.</p>
3	Dec 13, 2010	<p>Added Conformance Statement Overview and updated Introduction based on the new DICOM Conformance Statement template.</p> <p>Added reference to</p> <ul style="list-style-type: none"> • AW4.6 DICOM Conformance Statement • 22.0 Product Family DICOM Conformance Statement <p>Added support for</p> <ul style="list-style-type: none"> • X-Ray Radiation Dose SR <p>Updated Completion Flag and Verification Flag usage.</p> <p>Added changes after review:</p> <ul style="list-style-type: none"> • Added Issuer of Patient ID and Other Patient IDs Sequence in Patient Modules • Referenced Performed Procedure Step Sequence is not copied. • Verification flag is set to UNVERIFIED and Verifying Observer Sequence is cleared when a new report is saved. • Changed Mapping Resource to “99GEMS” for questionnaires in SR Document Content Module. • Series Description is copied for DICOM Encapsulated PDF in Encapsulated Document Series Module.



CONFORMANCE STATEMENT OVERVIEW

Reporting Tool is an Advantage Workstation application to display different SOP Classes of DICOM Structured Reports in various document formats (e.g. HTML, PDF, XML) with customizable layouts, and export them. Reporting Tool supports simple report editing where changes are saved as a new DICOM SR instance. Reporting Tool can create DICOM Secondary Capture Images from the PDF rendering of a report for DICOM archives that don't have DICOM SR support.

Table 0.1 provides an overview of the Storage SOP Classes supported by **Reporting Tool**.

Table 0.1 – SUPPORTED SOP CLASSES

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Object read / write		
Secondary Capture Image Storage	Yes	Yes
Basic Text SR	Yes	Yes
Enhanced SR	Yes	Yes
Comprehensive SR	Yes	Yes
Mammography CAD SR	Yes	Yes
Key Object Selection Document	Yes	Yes
X-Ray Radiation Dose SR	Yes	Yes
Encapsulated PDF Storage	Yes	Yes



TABLE OF CONTENTS

- 1. INTRODUCTION 6**
 - 1.1 OVERVIEW 6
 - 1.2 OVERALL DICOM CONFORMANCE STATEMENT DOCUMENT STRUCTURE 7
 - 1.3 INTENDED AUDIENCE 8
 - 1.4 SCOPE AND FIELD OF APPLICATION 8
 - 1.5 IMPORTANT REMARKS 8
 - 1.6 REFERENCES 9
 - 1.6.1 Platform Conformance Statements 9
 - 1.7 DEFINITIONS 9
- 2. NETWORK CONFORMANCE STATEMENT 11**
 - 2.1 INTRODUCTION 11
 - 2.2 IMPLEMENTATION MODEL 11
 - 2.3 AE SPECIFICATIONS 11
 - 2.4 COMMUNICATION PROFILES 12
 - 2.5 SUPPORT OF EXTENDED CHARACTER SETS 12
 - 2.6 CODES AND CONTROLLED TERMINOLOGY 12
 - 2.6.1 Standard coding scheme designators 12
 - 2.6.2 Private coding scheme designators 12
 - 2.7 SECURITY PROFILES 13
- 3. MEDIA STORAGE CONFORMANCE STATEMENT 14**
 - 3.1 INTRODUCTION 14
- 4. SECONDARY CAPTURE INFORMATION OBJECT IMPLEMENTATION 15**
 - 4.1 IOD MODULE TABLE 15
 - 4.2 INFORMATION MODULE DEFINITIONS 15
 - 4.2.1 Patient Module 16
 - 4.2.2 General Study Module 16
 - 4.2.3 General Series Module 16
 - 4.2.4 General Equipment Module 17
 - 4.2.5 SC Equipment Module 17
 - 4.2.6 General Image Module 17
 - 4.2.7 Image Pixel Module 17
 - 4.2.8 SC Image Module 18
 - 4.2.9 SOP Common Module 18
 - 4.3 Standard Extended and Private Data Attributes 18
 - 4.4 Standard Extended and Private CONTEXT Groups 18
- 5. BASIC TEXT, ENHANCED, COMPREHENSIVE, MAMMOGRAPHY CAD, X-RAY RADIATION DOSE SR AND KEY OBJECT SELECTION DOCUMENT INFORMATION OBJECT IMPLEMENTATION 19**
 - 5.1 IOD MODULE TABLE 19
 - 5.2 INFORMATION MODULE DEFINITIONS 20
 - 5.2.1 Patient Module 21
 - 5.2.2 General Study Module 21
 - 5.2.3 Patient Study Module 21
 - 5.2.4 SR Document Series Module 22
 - 5.2.5 General Equipment Module 22
 - 5.2.6 SR Document General Module 22
 - 5.2.7 SR Document Content Module 24
 - 5.2.8 SOP Common Module 26
 - 5.2.9 Key Object Document Series Module 27
 - 5.2.10 Key Object Document Module 27
 - 5.2.11 Synchronization Module 28
 - 5.2.12 Enhanced General Equipment Module 28
 - 5.3 Standard Extended and Private Data Attributes 28
 - 5.3.1 Private Group GEMS_0039 28
 - 5.4 Standard Extended and Private CONTEXT Groups 29
 - 5.4.1 Private Context Groups 29
 - 5.5 Standard Extended and Private Templates 30
 - 5.5.1 Private Templates 30



6.	ENCAPSULATED PDF INFORMATION OBJECT IMPLEMENTATION.....	32
6.1	IOD MODULE TABLE	32
6.2	INFORMATION MODULE DEFINITIONS	32
6.2.1	Patient Module	33
6.2.2	General Study Module	33
6.2.3	Encapsulated Document Series Module.....	33
6.2.4	General Equipment Module	34
6.2.5	SC Equipment Module.....	34
6.2.6	Encapsulated Document Module.....	34
6.2.7	SOP Common Module.....	34



1. INTRODUCTION

1.1 OVERVIEW

Reporting Tool software can display DICOM reports with customizable layouts and export them in different formats to various targets. It handles DICOM Structured Reports (SR) and DICOM Encapsulated PDFs. The tool is designed to run on:

- Advantage Workstation 4.2P, 4.3, 4.4, 4.5 and 4.6 (a Networked Medical Imaging Console).
- Signa HDXt and Discovery MR750 (MR Console).


This DICOM Conformance Statement is divided into Sections as described below:

- **Section 1 (Introduction)**, which describes the overall structure, intent, and references for this Conformance Statement
- **Section 2 (Network Conformance Statement)**, which specifies the GEHC equipment compliance to the DICOM requirements for the implementation of Networking features.
- **Section 3 (Media Storage Conformance Statement)**, which specifies the GEHC equipment compliance to the DICOM requirements for the implementation of Media Storage features.
- **Section 4 (Secondary Capture Information Object Implementation)**, which specifies the GEHC equipment compliance to DICOM requirements for the implementation of a Secondary Capture Information Object.
- **Section 5 (Structured Report Document Information Object Implementation)**, which specifies the GEHC equipment compliance to DICOM requirements for the implementation of a Structured Report Document Information Object.



1.2 OVERALL DICOM CONFORMANCE STATEMENT DOCUMENT STRUCTURE

The Documentation Structure of the GEHC DICOM Conformance Statements is shown in the Illustration below.



This document specifies the DICOM implementation. It is entitled:

Reporting Tool 2.7 (and up)
Conformance Statement for DICOM
Direction: **5215475-100**

This DICOM Conformance Statement documents the DICOM Conformance Statement and Technical Specification required to interoperate with the GEHC network interface.



The GEHC Conformance Statement, contained in this document, also specifies the Lower Layer communications, which it supports (e.g., TCP/IP). However, the Technical Specifications are defined in the DICOM Part 8 standard.

For more information regarding DICOM, copies of the Standard may be obtained on the Internet at <http://medical.nema.org>. Comments on the Standard may be addressed to:

DICOM Secretariat
NEMA
1300 N. 17th Street, Suite 1752
Rosslyn, VA 22209
USA
Phone : +1.703.841.3200

1.3 INTENDED AUDIENCE

The reader of this document is concerned with software design and/or system integration issues. It is assumed that the reader of this document is familiar with the DICOM Standard and with the terminology and concepts that are used in that Standard.

1.4 SCOPE AND FIELD OF APPLICATION

It is the intent of this document to provide an unambiguous specification for GEHC implementations. This specification, called a Conformance Statement, includes a DICOM Conformance Statement and is necessary to ensure proper processing and interpretation of GEHC medical data exchanged using DICOM. The GEHC Conformance Statements are available to the public.

The reader of this DICOM Conformance Statement should be aware that different GEHC devices are capable of using different Information Object Definitions. For example, a GEHC CT Scanner may send images using the CT Information Object, MR Information Object, Secondary Capture Object, etc.

Included in this DICOM Conformance Statement are the Module Definitions, which define all data elements, used by this GEHC implementation. If the user encounters unspecified private data elements while parsing a GEHC Data Set, the user is well advised to ignore those data elements (per the DICOM standard). Unspecified private data element information is subject to change without notice. If, however, the device is acting as a “full fidelity storage device”, it should retain and re-transmit all of the private data elements, which are sent by GEHC devices.

1.5 IMPORTANT REMARKS

The use of these DICOM Conformance Statements, in conjunction with the DICOM Standards, is intended to facilitate communication with GE imaging equipment. However, **by itself, it is not sufficient to ensure that inter-operation will be successful.** The **user (or user’s agent)** needs to proceed with caution and address at least four issues:

- **Integration** – The integration of any device into an overall system of interconnected devices goes beyond the scope of standards (DICOM v3.0), and of this introduction and associated DICOM Conformance Statements when interoperability with non-GE equipment is desired. The responsibility to analyze the applications requirements and to design a solution that integrates GE imaging equipment with non-GE systems is the **user’s** responsibility and should not be underestimated. The **user** is strongly advised to ensure that such an integration analysis is correctly performed.
- **Validation** – Testing the complete range of possible interactions between any GE device and non-GE devices, before the connection is declared operational, should not be overlooked. Therefore, the **user** should ensure that any non-GE provider accepts full responsibility for all validation required for their connection with GE devices. This includes the accuracy of the image data once it has crossed the interface between the GE imaging equipment and the non-GE device and the stability of the image data for the intended applications.
Such a validation is required before any clinical use (diagnosis and/or treatment) is performed. It applies when images acquired on GE imaging equipment are processed/displayed on a non-GE device, as well as when images acquired on non-GE equipment is processed/displayed on a GE console or workstation.
- **Future Evolution** – GE understands that the DICOM Standard will evolve to meet the user’s growing requirements. GE is actively involved in the development of the DICOM Standard. DICOM will incorporate new features and technologies and GE may follow the evolution of the Standard. The GEHC protocol is based on DICOM as specified in each DICOM Conformance Statement. Evolution of the Standard may require changes to devices, which have implemented DICOM. **In addition, GE reserves the right to discontinue or make changes to the support of communications features (on its**



products) described by these DICOM Conformance Statements. The user should ensure that any non-GE provider, which connects with GE devices, also plans for the future evolution of the DICOM Standard. Failure to do so will likely result in the loss of function and/or connectivity as the DICOM Standard changes and GE Products are enhanced to support these changes.

- **Interaction** – It is the sole responsibility of the non-GE provider to ensure that communication with the interfaced equipment does not cause degradation of GE imaging equipment performance and/or function.

1.6 REFERENCES

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>

1.6.1 Platform Conformance Statements

- Advantage Workstation 4.6 – 5404296-100
- Advantage Workstation 4.5 – 5324648-100
- Advantage Workstation 4.4 – 5181424-100
- Advantage Workstation 4.3 – 5138820-100
- Advantage Workstation 4.2P – 2381100-100
- Signa HDxt – DOC0489666
- Discovery MR750 – DOC0464555
- 22.0 Product Family DICOM Conformance Statement – DOC0700139

1.7 DEFINITIONS

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, and Computed Radiography Image Storage SOP Class.

Application Entity (AE) – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title – the externally known name of an *Application Entity*, used to identify a DICOM application to other DICOM applications on the network.

Application Context – the specification of the type of communication used between *Application Entities*. Example: DICOM network protocol.

Association – a network communication channel set up between *Application Entities*.

Attribute – a unit of information in an object definition; a data element identified by a *tag*. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004) and Procedure Code Sequence (0008,1032).

Information Object Definition (IOD) – the specified set of *Attributes* that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The *Attributes* may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG) – a set of standardized image compression techniques, available for use by DICOM applications.

Media Application Profile – the specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs)

Module – a set of *Attributes* within an *Information Object Definition* that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation – first phase of *Association* establishment that allows *Application Entities* to agree on the types of data to be exchanged and how that data will be encoded.

Presentation Context – the set of DICOM network services used over an *Association*, as negotiated between *Application Entities*; includes *Abstract Syntaxes* and *Transfer Syntaxes*.

Protocol Data Unit (PDU) – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.



Security Profile – a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an *Application Entity* to ensure confidentiality, integrity, and/or availability of exchanged DICOM data

Service Class Provider (SCP) – role of an *Application Entity* that provides a DICOM network service; typically, a server that performs operations requested by another *Application Entity* (*Service Class User*). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU) – role of an *Application Entity* that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

Service/Object Pair (SOP) Class – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

Service/Object Pair (SOP) Instance – an information object; a specific occurrence of information exchanged in a *SOP Class*. Examples: a specific x-ray image.

Tag – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

Transfer Syntax – the encoding used for exchange of DICOM information objects and messages. Examples: *JPEG* compressed (images), little endian explicit value representation.

Unique Identifier (UID) – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR) – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.



2. NETWORK CONFORMANCE STATEMENT

2.1 INTRODUCTION

This section of the DICOM Conformance Statement specifies the **Reporting Tool** compliance to DICOM requirements for **Networking** features.

Reporting Tool doesn't have own Networking features, these are provided by the underlying platform. For a complete description of the Networking features conformance refer to [Platform Conformance Statements](#).

2.2 IMPLEMENTATION MODEL

Refer to DICOM Conformance statement of the underlying platform.

2.3 AE SPECIFICATIONS

Reporting Tool creates new SR instances when:

- an existing report is modified and saved
- a patient questionnaire is filled out and saved.

Reporting Tool creates Secondary Capture (SC) Image Storage instances from the pages of a PDF report generated from a SR.

Reporting Tool creates Encapsulated PDF Storage instances from the PDF rendering of a SR.

Reporting Tool can display and update Key Object Selections but it doesn't create one from scratch.

2.3.1.1 SOP Classes

This Application Entity provides Standard Conformance to the following SOP Class(es):

TABLE 2.3-1
SOP Classes

SOP Class Name	SOP Class UID	Role
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11	SCU / SCP
Enhanced SR	1.2.840.10008.5.1.4.1.1.88.22	SCU / SCP
Comprehensive SR	1.2.840.10008.5.1.4.1.1.88.33	SCU / SCP
Mammography CAD SR	1.2.840.10008.5.1.4.1.1.88.50	SCU / SCP
Key Object Selection	1.2.840.10008.5.1.4.1.1.88.59	SCU / SCP
X-Ray Radiation Dose SR	1.2.840.10008.5.1.4.1.1.88.67	SCU / SCP
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	SCU / SCP
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	SCU / SCP

2.3.1.2 Association Policies

2.3.1.2.1 Implementation Identifying Information

TABLE 2.3-2
DICOM Implementation Class and Version

Implementation Class UID	1.2.840.113619.6.238
Implementation Version Name	RPT_<software version>

2.3.1.3 Association Initiation Policy

2.3.1.3.1 Presentation Context Table

These Transfer syntaxes are supported by the **Reporting Tool** application. For transfer syntaxes used for network transfer refer to [Platform Conformance Statements](#).

TABLE 2.3-3
Presentation Contexts



Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
X-Ray Radiation Dose SR	1.2.840.10008.5.1.4.1.1.88.67	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

2.4 COMMUNICATION PROFILES

Refer to [Platform Conformance Statements](#).

2.5 SUPPORT OF EXTENDED CHARACTER SETS

Reporting Tool supports the ISO_IR 100 (Latin alphabet No. 1) character set.

2.6 CODES AND CONTROLLED TERMINOLOGY

2.6.1 Standard coding scheme designators

This implementation makes use of the following standard coding scheme designators:

- **DCM**: DICOM Controlled Terminology; See [REFERENCES](#).
- **SRT**: SNOMED-RT (Referenced Terminology)
- **SNM3**: SNOMED Version 3 (used for backward compatibility)
- **UCUM**: [Unified Code for Units of Measure](#)

2.6.2 Private coding scheme designators

This implementation makes use of the following private coding scheme designators:



- 99GEMS

2.7 SECURITY PROFILES

The product does not conform to any defined DICOM Security Profiles.

It is assumed that the product is used within a secured environment. It is assumed that a secured environment includes at a minimum:

1. Firewall or router protections to ensure that only approved external hosts have network access to the product.
2. Firewall or router protections to ensure that the product only has network access to approved external hosts and services.
3. Any communications with external hosts and services outside the locally secured environment use appropriate secure network channels (such as a Virtual Private Network (VPN))



3. MEDIA STORAGE CONFORMANCE STATEMENT

3.1 INTRODUCTION

This section of the DICOM conformance statement specifies the **Reporting Tool** compliance to DICOM requirements for **Media Interchange**. It details the DICOM Media Storage Application Profiles and roles, which are supported by this product.

Reporting Tool doesn't have own Media Interchange implementation it is provided by the underlying platform. For a complete description of the Media Interchange conformance refer to [Platform Conformance Statements](#).



4. SECONDARY CAPTURE INFORMATION OBJECT IMPLEMENTATION

Reporting Tool

- creates Secondary Capture (SC) Image Storage instances from the pages of a PDF report generated from a SR.
- displays Secondary Capture (SC) Images that are referenced from an SR reading only the Image Pixel Module.

4.1 IOD MODULE TABLE

The Secondary Capture Information Object Definition comprises the modules of the following table, plus Standard Extended and Private attributes.

**TABLE 4.1-1
SC IMAGE IOD MODULES**

Information Entity	Module Name	Usage	Reference
Patient	Patient	Used	4.2.1
	Clinical Trial Subject	Not used	N/A
Study	General Study	Used	4.2.2
	Patient Study	Not used	N/A
	Clinical Trial Study	Not used	N/A
Series	General Series	Used	4.2.3
	Clinical Trial Series	Not used	N/A
Equipment	General Equipment	Used	4.2.4
	SC Equipment	Used	4.2.5
Image	General Image	Used	4.2.6
	Image Pixel	Used	4.2.7
	Device	Not used	N/A
	SC Image	Used	4.2.8
	Overlay Plane	Not used	N/A
	Modality LUT	Not used	N/A
	VOI LUT	Not used	N/A
	SOP Common	Used	4.2.9

4.2 INFORMATION MODULE DEFINITIONS

Please refer to DICOM v3.0 Standard Part 3 (Information Object Definitions) for a description of each of the entities and modules contained within the SC Information Object.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported. Type 1 & Type 2 Attributes are also included for completeness and to define what values they may take and where these values are obtained. It should be noted that they are the same ones as defined in the DICOM v3.0 Standard Part 3 (Information Object Definitions). Elements not listed in these modules are not generated.



4.2.1 Patient Module

**TABLE 4.2-1
PATIENT MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	Copied from source.
Patient ID	(0010,0020)	2	Copied from source.
Issuer of Patient ID	(0010,0021)	3	Copied from source.
Patient's Birth Date	(0010,0030)	2	Copied from source.
Patient's Sex	(0010,0040)	2	Copied from source. Enumerated Values: M = male F = female O = other
Other Patient IDs Sequence	(0010,1002)	3	Copied from source.
> Patient ID	(0010,0020)	1	Copied from source.
> Issuer of Patient ID	(0010,0021)	3	Copied from source.
> Type of Patient ID	(0010,0022)	1	Copied from source.

4.2.2 General Study Module

**TABLE 4.2-2
GENERAL STUDY MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
Study Instance UID	(0020,000D)	1	Copied from source.
Study Date	(0008,0020)	2	Copied from source.
Study Time	(0008,0030)	2	Copied from source.
Referring Physician's Name	(0008,0090)	2	Copied from source.
Study ID	(0020,0010)	2	Copied from source.
Accession Number	(0008,0050)	2	Copied from source.
Study Description	(0008,1030)	3	Copied from source.

4.2.3 General Series Module

**TABLE 4.2-3
GENERAL SERIES MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Enumerated Value: OT = Other
Series Instance UID	(0020,000E)	1	Generated with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.
Series Number	(0020,0011)	2	Copied from source.
Series Description	(0008,103E)	3	Copied from source.



4.2.4 General Equipment Module

TABLE 4.2-4
GENERAL EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	"GE MEDICAL SYSTEMS"
Institution Name	(0008,0080)	3	Copied from source.
Manufacturer's Model Name	(0008,1090)	3	"Reporting Tool"
Software Versions	(0018,1020)	3	Current software version.

4.2.5 SC Equipment Module

TABLE 4.2-5
SC EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Use
Conversion Type	(0008,0064)	1	SYN = Synthetic Image
Modality	(0008,0060)	3	Enumerated Value: OT = Other
Secondary Capture Device Manufacturer	(0018,1016)	3	"GE MEDICAL SYSTEMS"

4.2.6 General Image Module

TABLE 4.2-6
GENERAL IMAGE MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Instance Number	(0020,0013)	2	Generated.
Patient Orientation	(0020,0020)	2C	Empty.
Content Date	(0008,0023)	2C	Current date of creation.
Content Time	(0008,0033)	2C	Current time of creation.
Image Type	(0008,0008)	3	"DERIVED\SECONDARY"
Burned In Annotation	(0028,0301)	3	NO

4.2.7 Image Pixel Module

TABLE 4.2-7
IMAGE PIXEL MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Samples per Pixel	(0028,0002)	1	1 – for MONOCHROME2 3 – for RGB
Photometric Interpretation	(0028,0004)	1	MONOCHROME2 – when Samples per Pixel (0028,0002) has value 1 RGB – when Samples per Pixel (0028,0002) has value 3
Rows	(0028,0010)	1	Depends on PDF paper size and conversion resolution. E.g. with 72 dpi resolution: <ul style="list-style-type: none"> • 842 for A4 (210x297 mm) • 792 for Letter (8.5x11 in)



Columns	(0028,0011)	1	Depends on PDF paper size and conversion resolution. E.g. with 72 dpi resolution: <ul style="list-style-type: none"> • 595 for A4 (210x297 mm) • 612 for Letter (8.5x11 in)
Bits Allocated	(0028,0100)	1	8
Bits Stored	(0028,0101)	1	8
High Bit	(0028,0102)	1	7
Pixel Representation	(0028,0103)	1	0000H
Pixel Data	(7FE0,0010)	1	Derived from BMP generated from a PDF page.
Planar Configuration	(0028,0006)	1C	0 if Samples per Pixel (0028,0002) has a value greater than 1. Not present otherwise.

4.2.8 SC Image Module

**TABLE 4.2-8
SC IMAGE MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Use
Date of Secondary Capture	(0018,1012)	3	Current date of creation.
Time of Secondary Capture	(0018,1014)	3	Current time of creation.

4.2.9 SOP Common Module

**TABLE 4.2-9
SOP COMMON MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
SOP Class UID	(0008,0016)	1	Enumerated Values: 1.2.840.10008.5.1.4.1.1.7
SOP Instance UID	(0008,0018)	1	Generated with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.
Specific Character Set	(0008,0005)	1C	Copied from source if it contains non- empty value.
Instance Number	(0020,0013)	3	Generated.

4.3 STANDARD EXTENDED AND PRIVATE DATA ATTRIBUTES

Not used.

4.4 STANDARD EXTENDED AND PRIVATE CONTEXT GROUPS

Not used.



5. BASIC TEXT, ENHANCED, COMPREHENSIVE, MAMMOGRAPHY CAD, X-RAY RADIATION DOSE SR AND KEY OBJECT SELECTION DOCUMENT INFORMATION OBJECT IMPLEMENTATION

Reporting Tool creates new SR instances when:

- an existing report (from the supported SR SOP Classes) is modified and saved
- a patient questionnaire is created from scratch and saved.

Reporting Tool allows basic editing of the content of content items (e.g. basic text fields, TEXT items), but the user can't delete/insert one content item from/into the SR.

Reporting Tool allows report editing only when the Completion Flag of the SR is **PARTIAL**. I.e. reports that have Completion Flag of **COMPLETE** are finalized, cannot be edited.

5.1 IOD MODULE TABLE

Table 5.1-1 specifies the modules of **Basic Text, Enhanced, Comprehensive and Mammography CAD Structured Report** Information Object Definitions. SR specific modules are described in Section 5.2. Standard Extended and Private attributes are described in Section 5.3. The contents of the SR Document Content are constrained by the supported template, as identified in Section 5.2.7.1.2. Standard Extended and Private templates are further described in Section 5.5.

**TABLE 5.1-1
BASIC TEXT, ENHANCED, COMPREHENSIVE AND MAMMOGRAPHY CAD SR IOD MODULES**

Information Entity	Module	Usage	Reference
Patient	Patient	Used	5.2.1
	Clinical Trial Subject	Not used	N/A
Study	General Study	Used	5.2.2
	Patient Study	Used	5.2.3
	Clinical Trial Study	Not used	N/A
Series	SR Document Series	Used	5.2.4
	Clinical Trial Series	Not used	N/A
Equipment	General Equipment	Used	5.2.5
Document	SR Document General	Used	5.2.6
	SR Document Content	Used	5.2.7
	SOP Common	Used	5.2.8

Table 5.1-2 specifies the modules of **Key Object Selection Document** Information Object Definition.



**TABLE 5.1-2
KEY OBJECT SELECTION IOD MODULES**

Information Entity	Module	Usage	Reference
Patient	Patient	Used	5.2.1
	Clinical Trial Subject	Not used	N/A
Study	General Study	Used	5.2.2
	Patient Study	Used	5.2.3
	Clinical Trial Study	Not used	N/A
Series	Key Object Document Series	Used	5.2.9
	Clinical Trial Series	Not used	N/A
Equipment	General Equipment	Used	5.2.5
Document	Key Object Document	Used	5.2.10
	SR Document Content	Used	5.2.7
	SOP Common	Used	5.2.8

Table 5.1-3 specifies the modules of **X-Ray Radiation Dose Structured Report** Information Object Definition.

**TABLE 5.1-3
X-RAY RADIATION DOSE SR IOD MODULES**

Information Entity	Module	Usage	Reference
Patient	Patient	Used	5.2.1
	Clinical Trial Subject	Not used	N/A
Study	General Study	Used	5.2.2
	Patient Study	Used	5.2.3
	Clinical Trial Study	Not used	N/A
Series	SR Document Series	Used	5.2.4
	Clinical Trial Series	Not used	N/A
Frame of Reference	Synchronization	Used if present	5.2.11
Equipment	General Equipment	Used	5.2.5
	Enhanced General Equipment	Used	5.2.12
Document	SR Document General	Used	5.2.6
	SR Document Content	Used	5.2.7
	SOP Common	Used	5.2.8

5.2 INFORMATION MODULE DEFINITIONS

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the SR Information Objects.

Displaying of data elements is configurable i.e. they can be rendered in the report or stay hidden. The “Displayed” column in the following tables describe the elements Reporting Tool can display. Elements not listed in these modules are not used.



5.2.1 Patient Module

TABLE 5.2-1
PATIENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Patient's Name	(0010,0010)	2	Copied from source.	Yes
Patient ID	(0010,0020)	2	Copied from source.	Yes
Issuer of Patient ID	(0010,0021)	3	Copied from source.	
Patient's Birth Date	(0010,0030)	2	Copied from source.	Yes
Patient's Sex	(0010,0040)	2	Copied from source. Enumerated Values: M = male F = female O = other	Yes
Patient's Birth Time	(0010,0032)	3	Copied from source.	No
Other Patient IDs Sequence	(0010,1002)	3	Copied from source.	No
> Patient ID	(0010,0020)	1	Copied from source.	No
> Issuer of Patient ID	(0010,0021)	3	Copied from source.	No
> Type of Patient ID	(0010,0022)	1	Copied from source.	No
Ethnic Group	(0010,2160)	3	Copied from source or entered by the user.	Yes
Patient Comments	(0010,4000)	3	Copied from source.	No

5.2.2 General Study Module

TABLE 5.2-2
GENERAL STUDY MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Study Instance UID	(0020,000D)	1	Copied from source.	No
Study Date	(0008,0020)	2	Copied from source.	Yes
Study Time	(0008,0030)	2	Copied from source.	Yes
Referring Physician's Name	(0008,0090)	2	Copied from source or entered by the user.	Yes
Study ID	(0020,0010)	2	Copied from source.	Yes
Accession Number	(0008,0050)	2	Copied from source or entered by the user.	Yes
Study Description	(0008,1030)	3	Copied from source or empty.	No

5.2.3 Patient Study Module

TABLE 5.2-3
PATIENT STUDY MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Patient's Age	(0010,1010)	3	Copied from source or empty.	Yes
Patient's Size	(0010,1020)	3	Copied from source or entered by the user.	Yes
Patient's Weight	(0010,1030)	3	Copied from source or entered by the user.	Yes
Occupation	(0010,2180)	3	Copied from source or empty.	No
Additional Patient's History	(0010,21B0)	3	Copied from source or empty.	No



5.2.4 SR Document Series Module

TABLE 5.2-4
SR DOCUMENT SERIES MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Modality	(0008,0060)	1	Copied from source. Enumerated Value: SR = SR Document	No
Series Instance UID	(0020,000E)	1	Copied from source or generated if the SR is created from scratch with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.	No
Series Number	(0020,0011)	1	Copied from source or set to "1" if the SR is created from scratch.	No
Series Description	(0008,103E)	3	Copied from source.	No
Referenced Performed Procedure Step Sequence	(0008,1111)	2	Empty.	No

5.2.5 General Equipment Module

TABLE 5.2-5
GENERAL EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Manufacturer	(0008,0070)	2	Copied from source or "GE MEDICAL SYSTEMS" if the SR is created from scratch.	No
Institution Name	(0008,0080)	3	Copied from source or empty.	No
Institution Address	(0008,0081)	3	Copied from source or empty.	No
Station Name	(0008,1010)	3	Copied from source or empty.	No
Institutional Department Name	(0008,1040)	3	Copied from source or empty.	No
Manufacturer's Model Name	(0008,1090)	3	Copied from source or "Reporting Tool".	No
Device Serial Number	(0018,1000)	3	Copied from source or empty.	No

5.2.6 SR Document General Module

TABLE 5.2-6
SR DOCUMENT GENERAL MODULE ATTRIBUTES



Attribute Name	Tag	Type	Attribute Description	Displayed
Instance Number	(0020,0013)	1	Copied from source or set to "1" if the SR is created from scratch.	No
Completion Flag	(0040,A491)	1	Copied from source, selected by the user or PARTIAL if the SR is created from scratch. Enumerated Values: PARTIAL = Partial content. COMPLETE = Complete content.	Yes
Completion Flag Description	(0040,A492)	3	Copied from source or empty.	No
Verification Flag	(0040,A493)	1	Set to UNVERIFIED.	Yes
Content Date	(0008,0023)	1	Current date of creation.	Yes
Content Time	(0008,0033)	1	Current time of creation.	Yes
Verifying Observer Sequence	(0040,A073)	1C	Not copied.	Yes
>Verifying Observer Name	(0040,A075)	1	Not copied.	Yes
>Verifying Observer Identification Code Sequence	(0040,A088)	2	Not copied.	No
>Verifying Organization	(0040,A027)	1	Not copied.	Yes
>Verification DateTime	(0040,A030)	1	Not copied.	Yes
Predecessor Documents Sequence	(0040,A360)	1C	Reference to source SR when an existing report is amended. Not present if the SR is created from scratch.	Yes
>Study Instance UID	(0020,000D)	1	The Study Instance UID of the source SR.	No
>Referenced Series Sequence	(0008,1115)	1		No
>>Series Instance UID	(0020,000E)	1	The Series Instance UID of the source SR.	No
>>Retrieve AE Title	(0008,0054)	3	Empty.	No
>>Storage Media File-Set ID	(0088,0130)	3	Empty.	No
>>Storage Media File-Set UID	(0088,0140)	3	Empty.	No
>>Referenced SOP Sequence	(0008,1199)	1		No
>>>Referenced SOP Class UID	(0008,1150)	1	The SOP Class UID of the source SR.	No
>>>Referenced SOP Instance UID	(0008,1155)	1	The SOP Instance UID of the source SR	No
Identical Documents Sequence	(0040,A525)	1C	Not present.	No
Referenced Request Sequence	(0040,A370)	1C	Not present.	No
Performed Procedure Code Sequence	(0040,A372)	2	Copied from source or empty if the SR is created from scratch.	No
Current Requested Procedure Evidence Sequence	(0040,A375)	1C	Copied from source or not present if the SR is created from scratch.	No



Pertinent Other Evidence Sequence	(0040,A385)	1C	Copied from source or not present if the SR is created from scratch.	No
-----------------------------------	-------------	----	--	----

5.2.7 SR Document Content Module

TABLE 5.2-7
SR DOCUMENT CONTENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
<i>Include 'Document Content Macro' Table 5.2-8 with a Value Type (0040,A040) of CONTAINER.</i>				
<i>Include 'Document Relationship Macro' Table 5.2-9.</i>				

TABLE 5.2-8
DOCUMENT CONTENT MACRO ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Value Type	(0040,A040)	1	Copied from source or CONTAINER for the root content item if the SR is created from scratch. Defined Terms: TEXT NUM CODE DATETIME DATE TIME UIDREF PNAME COMPOSITE IMAGE WAVEFORM SCoord TCoord CONTAINER	No
Concept Name Code Sequence	(0040,A043)	1C	Conveys Document Title. Copied from source or set to the triplet below for patient questionnaires created from scratch.	Yes
> Code Value	(0008,0100)	1C	Copied from source or "PQ-100" for patient questionnaires created from scratch.	No
> Coding Scheme Designator	(0008,0102)	1C	Copied from source or "99GEMS" for patient questionnaires created from scratch.	No
> Code Meaning	(0008,0104)	1C	Copied from source or "Questionnaire" for patient questionnaires created from scratch.	Yes
<i>Include 'Container Macro' Table 5.2-10 if and only if Value Type (0040,A040) is CONTAINER.</i>				

TABLE 5.2-9
DOCUMENT RELATIONSHIP MACRO ATTRIBUTES



Attribute Name	Tag	Type	Attribute Description	Displayed
Observation DateTime	(0040,A032)	1C	Copied from source or current date and time.	No
Content Sequence	(0040,A730)	1C	SR content, the sequence of top-level content items.	Yes
> Relationship Type	(0040,A010)	1	Copied from source. Defined Terms: CONTAINS HAS PROPERTIES HAS OBS CONTEXT HAS ACQ CONTEXT INFERRED FROM SELECTED FROM HAS CONCEPT MOD	No
> Referenced Content Item Identifier	(0040,DB73)	1C	Not used. Note: Comprehensive SRs are displayed and editable but the references between the content items are not handled.	No

**TABLE 5.2-10
CONTAINER MACRO ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description	Displayed
Continuity of Content	(0040,A050)	1	Copied from source or "SEPARATE" if the SR is created from scratch. Enumerated Values: SEPARATE CONTINUOUS	No
Content Template Sequence	(0040,A504)	1C	Copied from source or filled in for Questionnaire.	No
>Mapping Resource	(0008,0105)	1	"99GEMS" for Questionnaire.	No
>Template Identifier	(0040,DB00)	1	"PQ_100" for Questionnaire.	No

5.2.7.1 SR Document Content Descriptions

5.2.7.1.1 Value Type rendering and editing

Table 5.2-11 describes how the different content item types are rendered in a generic format and which one is editable. The generic layout is an indented display of the nested content items' recursive traversal. Each content item's **Code Meaning** (0008,0104) in the corresponding Concept Name Code Sequence (0040,A043) is displayed together with the actual value itself. Reference Value Types are not editable, which means that it's not possible to add, update or remove reference to an IMAGE, WAVEFORM, COMPOSITE or UIDREF content item.

**TABLE 5.2-11
RENDERING AND EDITING BY VALUE TYPES**



Value Type (0040,A040)	Rendering	Editable
TEXT	Text Value (0040,A160)	Yes
NUM	Numeric Value (0040,A30A) and Code Value (0008,0100) from Measurement Units Code Sequence (0040,08EA)	No
CODE	Code Meaning (0008,0104) of Concept Code Sequence (0040,A168)	No
DATETIME	Date Time (0040,A120) in YYYY-MM-DD, hh:mm:ss format	No
DATE	Date (0040,A121) in YYYY-MM-DD format	No
TIME	Time (0040,A122) in hh:mm:ss format	No
UIDREF	UID (0040,A124)	No
PNAME	Person Name (0040,A123) in name_prefix given_name_complex middle_name family_name_complex name_suffix format	No
COMPOSITE	Referenced SOP Instance UID (0008,1155)	No
IMAGE	JPEG image of the DICOM object identified by Referenced SOP Instance UID (0008,1155). GSPS is not supported. Multi-frame image is not supported.	No
WAVEFORM	Referenced SOP Instance UID (0008,1155)	No
SCCOORD	Not rendered	N/A
TCOORD	Not rendered	N/A
CONTAINER	Recursive rendering of child content items	N/A

5.2.7.1.2 Content Template

The product supports the following root Templates.

**TABLE 5.2-12
ROOT TEMPLATES**

SOP Class	Template ID	Template Name	Use
Basic Text SR	Any		Display / Update
Enhanced SR	Any		Display / Update
Comprehensive SR	Any		Display / Update
Mammography CAD SR	TID 4000	Mammography CAD Document Root	Display / Update
Key Object Selection Document	TID 2010	Key Object Selection	Display / Update
X-Ray Radiation Dose SR	TID 10001	Projection X-Ray Radiation Dose	Display / Update
	TID 10011	CT Radiation Dose	Display / Update

5.2.8 SOP Common Module

TABLE 5.2-13



SOP COMMON MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
SOP Class UID	(0008,0016)	1	Enumerated Values: 1.2.840.10008.5.1.4.1.1.88.11 1.2.840.10008.5.1.4.1.1.88.22 1.2.840.10008.5.1.4.1.1.88.33 1.2.840.10008.5.1.4.1.1.88.50 1.2.840.10008.5.1.4.1.1.88.59 1.2.840.10008.5.1.4.1.1.88.67
SOP Instance UID	(0008,0018)	1	Generated with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.
Specific Character Set	(0008,0005)	1C	Copied from source or "ISO_IR 100".
Instance Creation Date	(0008,0012)	3	Current date of creation.
Instance Creation Time	(0008,0013)	3	Current time of creation.
Instance Creator UID	(0008,0014)	3	Empty.

5.2.9 Key Object Document Series Module

**TABLE 5.2-14
KEY OBJECT DOCUMENT SERIES MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description	Displayed
Modality	(0008,0060)	1	Copied from source. Enumerated Value: KO = Key Object Selection	No
Series Instance UID	(0020,000E)	1	Copied from source.	No
Series Number	(0020,0011)	1	Copied from source.	No
Series Description	(0008,103E)	3	Copied from source.	No
Referenced Performed Procedure Step Sequence	(0008,1111)	2	Empty.	No

5.2.10 Key Object Document Module

**TABLE 5.2-15
KEY OBJECT DOCUMENT MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description	Displayed
Instance Number	(0020,0013)	1	Copied from source.	No
Content Date	(0008,0023)	1	Current date of creation.	Yes
Content Time	(0008,0033)	1	Current time of creation.	Yes
Referenced Request Sequence	(0040,A370)	1C	Not present.	No
Current Requested Procedure Evidence Sequence	(0040,A375)	1	Copied from source.	No
Identical Documents Sequence	(0040,A525)	1C	Not present.	No



5.2.11 Synchronization Module

TABLE 5.2-16
SYNCHRONIZATION MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Synchronization Frame of Reference UID	(0020,0200)	1	Copied from source.	No
Synchronization Trigger	(0018,106A)	1	Copied from source.	No
Trigger Source or Type	(0018,1061)	3	Copied from source.	No
Synchronization Channel	(0018,106C)	1C	Copied from source.	No
Acquisition Time Synchronized	(0018,1800)	1	Copied from source.	No
Time Source	(0018,1801)	3	Copied from source.	No
Time Distribution Protocol	(0018,1802)	3	Copied from source.	No
NTP Source Address	(0018,1803)	3	Copied from source.	No

5.2.12 Enhanced General Equipment Module

TABLE 5.2-17
ENHANCED GENERAL EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description	Displayed
Manufacturer	(0008,0070)	1	Copied from source.	No
Manufacturer's Model Name	(0008,1090)	1	Copied from source or "Reporting Tool".	No
Device Serial Number	(0018,1000)	1	Copied from source.	No
Software Versions	(0018,1020)	1	Copied from source.	No

5.3 STANDARD EXTENDED AND PRIVATE DATA ATTRIBUTES

The Product supports the Standard and Private Attributes defined in the following sections in Standard Extended SR SOP Instances as Type 3 data elements.

5.3.1 Private Group GEMS_0039

TABLE 5.3-1
PRIVATE GROUP GEMS_0039

Attribute Name	Tag	VR	VM	Attribute Description and Use
Private Creator Identification	(0039,0010)	LO	1	"REPORT_FROM_APP"
Application specific data	(0039,1095)	LO	1	VV#<application_version>#<application_name>

Reporting Tool is expecting only (0039,1095).

This data is used to render the report in application specific format, which is different from the generic format described in 5.2.7.1.1.



5.4 STANDARD EXTENDED AND PRIVATE CONTEXT GROUPS

5.4.1 Private Context Groups

5.4.1.1 CID PQ-010 Lung Analysis Questions

Context ID PQ-010

Lung Analysis Questions

Type: Extensible

Version: 20070925

Coding Scheme Designator (0008,0102)	Coding Scheme Version (0008,0103)	Code Value (0008,0100)	Code Meaning (0008,0104)
99GEMS		PQ-1000	Reason for exam
99GEMS		PQ-1010	Pre-exam comments
99GEMS		PQ-1020	Currently smoking cigarettes
99GEMS		PQ-1030	If not currently smoking, when did patient quit?
99GEMS		PQ-1040	Age when first started smoking cigarettes
99GEMS		PQ-1050	How many times has patient tried to quit smoking?
99GEMS		PQ-1060	How many pack(s) a day has the patient smoked?
99GEMS		PQ-1070	For how many years has the patient smoked?
99GEMS		PQ-1080	Number of cigarette pack years
99GEMS		PQ-1090	Select other tobacco products used
99GEMS		PQ-1100	Has patient ever used smokeless tobacco products?
99GEMS		PQ-1105	If so, for how many years?
99GEMS		PQ-1110	Has patient been exposed to significant second-hand smoke in his/her home?
99GEMS		PQ-1120	Has patient been exposed to significant second-hand smoke in his/her workplace?
99GEMS		PQ-1130	Has patient been exposed to significant second-hand smoke in social settings?
99GEMS		PQ-1140	Do/did patient's parents smoke?
99GEMS		PQ-1150	Does patient's spouse smoke?
99GEMS		PQ-1160	Has patient ever worked with asbestos?
99GEMS		PQ-1170	Has patient ever worked in a shipyard or with pipe insulation?
99GEMS		PQ-1180	Has patient ever worked with heavy metals such as lead, arsenic, or mercury?
99GEMS		PQ-1190	Has patient ever worked in a mine?
99GEMS		PQ-1200	Has patient had significant exposure to radon gas?
99GEMS		PQ-1210	Does the patient have a past history of cancer?
99GEMS		PQ-1220	Respiratory symptoms
99GEMS		PQ-1230	Has the patient ever been told that he/she had an abnormal chest x-ray?
99GEMS		PQ-1240	If so, what was the diagnosis?
99GEMS		PQ-1250	Has patient ever had a positive skin test for tuberculosis?
99GEMS		PQ-1260	Does patient have a history of asthma?
99GEMS		PQ-1270	Does the patient have a history of chronic or frequent bronchitis?
99GEMS		PQ-1280	Does the patient have a chronic smoker's cough?
99GEMS		PQ-1290	Has the patient had a fungal infection of the lungs?
99GEMS		PQ-1300	Has the patient ever been diagnosed with emphysema?
99GEMS		PQ-1310	Family history of lung cancer
99GEMS		PQ-1320	Please check any conditions for which the patient has seen a health care provider or taken medication in the past five years



Coding Scheme Designator (0008,0102)	Coding Scheme Version (0008,0103)	Code Value (0008,0100)	Code Meaning (0008,0104)
99GEMS		PQ-1330	Current living location
SRT		J-00000	Occupation

5.4.1.2 CID PQ-020 Oncology Questions

Context ID PQ-020

Oncology Questions

Type: Extensible

Version: 20070925

Coding Scheme Designator (0008,0102)	Coding Scheme Version (0008,0103)	Code Value (0008,0100)	Code Meaning (0008,0104)
99GEMS		PQ-2000	Cancer diagnosis
99GEMS		PQ-2010	AJCC/UICC Staging
SRT		G-F150	T category
SRT		R-40030	N category
SRT		R-40031	M category
99GEMS		PQ-2020	Stage grouping
99GEMS		PQ-2030	Medical and surgical history
99GEMS		PQ-2040	Current medical and surgical history
99GEMS		PQ-2050	Allergies
99GEMS		PQ-2060	Dermatologic status
99GEMS		PQ-2070	Laboratory data history
99GEMS		PQ-2080	Treatment history

5.5 STANDARD EXTENDED AND PRIVATE TEMPLATES

5.5.1 Private Templates

5.5.1.1 TID PQ_100 Questionnaire Document Root Template

This template forms the top of a content tree for questionnaire reports generated by Reporting Tool.

TID PQ_100

QUESTIONNAIRE DOCUMENT ROOT

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (PQ-100, 99GEMS, "Questionnaire")	1	M		
2	>	CONTAINS	INCLUDE	DTID (PQ_200) Patient Questions	1	M		

5.5.1.2 TID PQ_200 Patient Questions Template

This template describes the list of questions in the questionnaire.



TID PQ_200
PATIENT QUESTIONS

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (PQ-0110, 99GEMS, "Patient Questions")	1	M		
2	>	CONTAINS	INCLUDE	DTID (PQ_300) Patient Question Answer	1-n	M		

5.5.1.3 TID PQ_300 Patient Question Answer Template

This template describes one question with its corresponding answer.

TID PQ_300
PATIENT QUESTION ANSWER

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1		CONTAINS	CODE	EV (PQ-0120, 99GEMS, "Question")	1-n	MC	XOR row 4	DCID(PQ-010) DCID(PQ-020)
2	>	HAS CONCEPT MOD	TEXT	EV (121051,DCM,"Equivalent Meaning of Value")	1	M		
3	>	HAS PROPERTIES	TEXT	EV (PQ-0130, 99GEMS, "Answer")	1	U		
4		CONTAINS	TEXT	EV (PQ-0120, 99GEMS, "Question")	1-n	MC	XOR row 1	
5	>	HAS PROPERTIES	TEXT	EV (PQ-0130, 99GEMS, "Answer")	1	U		



6. ENCAPSULATED PDF INFORMATION OBJECT IMPLEMENTATION

Reporting Tool creates Encapsulated PDF Storage instances from the PDF rendering of a SR (source).

6.1 IOD MODULE TABLE

Table 6.1-1 specifies the modules of **Encapsulated PDF** Information Object Definition.

TABLE 6.1-1
ENCAPSULATED PDF IOD MODULES

Information Entity	Module	Usage	Reference
Patient	Patient	Used	6.2.1
	Clinical Trial Subject	Not used	N/A
Study	General Study	Used	6.2.2
	Patient Study	Not used	N/A
	Clinical Trial Study	Not used	N/A
Series	Encapsulated Document Series	Used	6.2.3
	Clinical Trial Series	Not used	N/A
Equipment	General Equipment	Used	6.2.4
	SC Equipment	Used	6.2.5
Encapsulated Document	Encapsulated Document	Used	6.2.6
	SOP Common	Used	6.2.7

6.2 INFORMATION MODULE DEFINITIONS

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the Encapsulated PDF Information Object.



6.2.1 Patient Module

TABLE 6.2-1
PATIENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	Copied from source.
Patient ID	(0010,0020)	2	Copied from source.
Issuer of Patient ID	(0010,0021)	3	Copied from source.
Patient's Birth Date	(0010,0030)	2	Copied from source.
Patient's Sex	(0010,0040)	2	Copied from source. Enumerated Values: M = male F = female O = other
Other Patient IDs Sequence	(0010,1002)	3	Copied from source.
> Patient ID	(0010,0020)	1	Copied from source.
> Issuer of Patient ID	(0010,0021)	3	Copied from source.
> Type of Patient ID	(0010,0022)	1	Copied from source.

6.2.2 General Study Module

TABLE 6.2-2
GENERAL STUDY MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Study Instance UID	(0020,000D)	1	Copied from source.
Study Date	(0008,0020)	2	Copied from source.
Study Time	(0008,0030)	2	Copied from source.
Referring Physician's Name	(0008,0090)	2	Copied from source.
Study ID	(0020,0010)	2	Copied from source.
Accession Number	(0008,0050)	2	Copied from source.

6.2.3 Encapsulated Document Series Module

TABLE 6.2-3
ENCAPSULATED DOCUMENT SERIES MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Enumerated Value: OT = Other
Series Instance UID	(0020,000E)	1	Generated with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.
Series Number	(0020,0011)	1	Copied from source.
Series Description	(0008,103E)	3	Copied from source.



6.2.4 General Equipment Module

TABLE 6.2-4
GENERAL EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	"GE MEDICAL SYSTEMS"

6.2.5 SC Equipment Module

TABLE 6.2-5
SC EQUIPMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Use
Conversion Type	(0008,0064)	1	WSD = Workstation

6.2.6 Encapsulated Document Module

TABLE 6.2-6
ENCAPSULATED DOCUMENT MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Instance Number	(0020,0013)	1	Generated.
Content Date	(0008,0023)	2	Copied from source.
Content Time	(0008,0033)	2	Copied from source.
Acquisition DateTime	(0008,002A)	2	Copied from source.
Burned In Annotation	(0028,0301)	1	Copied from source if it is Encapsulated PDF. If the source is SR and the PDF rendering contains sufficient burned in annotation to identify the patient and date the data was acquired then YES, otherwise NO.
Source Instance Sequence	(0042,0013)	1C	Copied from source or filled with SOP Instance Reference from the source.
Document Title	(0042,0010)	2	Copied from source or filled with Code Meaning from source Concept Name Code Sequence.
Concept Name Code Sequence	(0040,A043)	2	Copied from source.
MIME Type of Encapsulated Document	(0042,0012)	1	"application/pdf"
Encapsulated Document	(0042,0011)	1	Byte stream of PDF rendering of source.

6.2.7 SOP Common Module

TABLE 6.2-7
SOP COMMON MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
SOP Class UID	(0008,0016)	1	Enumerated Values: 1.2.840.10008.5.1.4.1.1.104.1
SOP Instance UID	(0008,0018)	1	Generated with format 1.2.840.113619.2.238.id where id is a unique identifier of the instance with station information and timestamp.
Specific Character Set	(0008,0005)	1C	Copied from source or "ISO_IR 100".