

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT (Equipment locations, heat loads, component weights, environmental specs)	A1
STRUCTURAL LAYOUT (Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	S1
STRUCTURAL DETAILS (Floor and Ceiling loading information)	S2
ELECTRICAL LAYOUT (Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	E1
ELECTRICAL SPECIFICATIONS (Maximum wiring run lengths, interconnect diagram, system power specifications)	E2
ELECTRICAL DETAILS	E3
EQUIPMENT DETAILS	D1 THRU D2

These equipment IS drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the IS and operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

* REQUIRED REFERENCE *

LightSpeed RT 16 / Xtra
Pre Installation Manual
5177460-100

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



CT Site Planning



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist						
GEHC Global Order #:			Customer:			
GEHC On-site Representative:			MI Supplier:			
Name of customer reviewed with:			Lead Installer:			
GEHC PMI#:			Phone Number:			
Target Site Prep Completion Date:			Helper:			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.						
For MR Magnet Delivery: Ensure cryogen vents, power for the cooling system and exhaust fan system are installed and operational (0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20) and chilled water supply is available 24x7 that meets system cooling equipment requirements. Broadband/phone line available for magnet monitor.						
Item #	GEHC Minimum Requirements	Storage: Is item ready?	Protect (Pre ship): Is this item ready? Will item be ready?	Verify (Delivery): Is item ready?	Validate (Mech Install): Is item ready?	Comments If "N", please enter in comments or action plan
1	Equipment installation drawings must show clearances, room size, equipment placement and must meet clearance requirements. Deviations that meet installation requirements may be red-lined, if allowed by local code. Seismic requirements identified on construction drawings. In the following states-NC, SC, AR, DC, WA - Verify State approved shielding plans match GE drawing room dimensions and requirements.					
2	Delivery route to installation or storage area meets requirements and has been discussed and scheduled with the customer. Ensure floor protection is discussed, requirements identified, and will be available at time of delivery and installation.					
3	Rooms that will contain equipment including storage areas-not in scan suite, are dust free. Provisions taken to maintain a dust free room. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility.					
4	In room HVAC ductwork and units (in room) must be mechanically installed and dust free. Installation rooms appear to meet environmental conditions (see Further Definitions) and observed issues have been communicated to the customer. If being stored, storage area-not in scan suite, must meet PM storage criteria.					
5	Ceiling grid is installed. Permanent lighting is installed and operational. Unistrut (or equivalent) location and spacing was measured and is consistent with the requirements of the installation drawings.					
6	Floor is clean and prepared for final floor covering. For MR, CT & Nuc scan rooms, floor levelness was measured and does not exceed tolerances specified in GEHC's applicable PM, and no visible floor surface defects were observed.					
7	Access to a working phone at the facility for emergency use, including MR magnet delivery.					
8	All walls primed (final coat not needed on Day 1)					
9	Mechanical supplier has been provided with a set of equipment installation drawings for reference. For California, permitted construction drawings or PMI-specified installation drawings are required.					
10	Conduit/electrical cable ducting/dividers/ access flooring installed, with the exception of surface-mounted floor ducting. Wiring to the main disconnect panel is installed and compliant with equipment installation drawings or pre-installation manual.					

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin
Copyright 2009 General Electric Company - Proprietary to GE

SHEET TITLE: SITE READINESS
MODALITY TYPE: LIGHTSPEED RT16/XTRA
THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IT IS NOT TO BE USED FOR CONSTRUCTION PURPOSES. THE CUSTOMER IS RESPONSIBLE FOR OBTAINING ALL NECESSARY PERMITS AND APPROVALS FOR ACTUAL CONSTRUCTION PURPOSES. HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
6-64f
TYPICAL FINAL
(with GT1700 Table)

PROJECT	REVISION
6-64f	05
DATE:	25 Jan. 12
DRAWN BY:	JGA
CHECKED BY:	JGA

REVISION HISTORY:

SHEET
C1

PIM RT1

GE EQUIPMENT LISTING

EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER : NEITHER A QUOTE OR CON WAS ISSUED AT THE DATE OF THESE DRAWINGS

NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.

ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	EQUIPMENT CROSS REFERENCE CHART		STRC PLAN	ELEC PLAN	
							SSMIC STATUS	P = PREAPPROVAL C = CALCULATIONS/ PENDING APPROVAL S = SPECIFICATIONS ONLY			
1	1		INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	22 lbs		B5031	B5031S	IC	IC		S
2	1		INJECTOR CONTROL CONSOLE	8 lbs		B8007C		IC	IC		S
3	1		MAIN UNIT	13 lbs		B8007C		IE	IE		-
4	1		UPS SYSTEM	350 lbs		B7999ZA		UPS			-
6	1		POWER DISTRIBUTION UNIT	701 lbs	3402 btu	B7858D		PM	PM		S
7	1		OPERATOR'S CONSOLE / COMPUTER	493 lbs	8126 btu	B7858A		OC	OC		S
8	1		OPERATOR'S CHAIR								-
9	1		LCD MONITOR	61 lbs	300 btu	B5031S			SVM		S
10	1		STORAGE CABINET (EMPTY CABINET WEIGHT)	90 lbs		M3300S					-
11	1		REAR CABLE COVER								-
11	1		LIGHTSPEED RT16/XTRA	3891 lbs	43010 btu	B7817B B7817A B7817C B7817D B7800C B7864B		B7817	CTT		C
12	1		GT1700 PATIENT TABLE WITH EXTENDED TABLE TOP.	881 lbs	703 btu	B7917					-
13	1		ADVANTAGE WORKSTATION WITH TWO LCD MONITORS	81 lbs	1109 btu	M1013AW					S

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.


50	2		PATIENT POSITIONING LASER LIGHT	66 lbs		B7998B	B7998D				S
51	1		VARIAN CAMERA HOLDER (WHEN CAMERA IS IN STANDBY POSITION)								-
52	1		VARIAN POWER SUPPLY MODULE (GE PROJ MGR. PREINSTALL ITEM. PRE-ORDER FROM VARIAN)	13 lbs		R4504D		PSM			-
53	1		VARIAN OPERATORS WORKSTATION	48 lbs	819 btu			RWS			-
54	1		VARIAN CAMERA AND BRACKET	8 lbs	341 btu	E8819KA		RGC			-
55	1		WORKSTATION COMPUTER					CLC			S
56	1		PATIENT POSITIONING LASER LIGHT	66 lbs		B7998A	B7998C	CL			S
57	1		KEYPAD	2 lbs		B7998F		KP			S
58	1		MAIN DISCONNECT CONTROL	132 lbs		E4502AE		A1			C

SCALE: 1/4" = 1'-0"

EQUIPMENT LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-0"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	WORKSTATION TABLE
61	COUNTER TOP WITH SINK, BASE AND WALL CABINETS
62	LEAD GLASS WINDOW
63	COUNTER TOP FOR EQUIPMENT - PROVIDE GROMMETTED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO REARWAY BELOW COUNTERTOP.
64	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL: 800-200-9760 GE CAT. NO. X1ABW-DF-XIU
65	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W x 83 IN. H (1118mm x 2108mm), CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT CONTROL PANEL REFERENCE JUNCTION POINT - WLS ON SHEET 611 FOR DETAILED DESCRIPTION - E4506RL FOR WARNING LIGHT CONTROL ONLY.
----	---

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC IS SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PRECIPITATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC.
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM

SITE ENVIRONMENT SPECIFICATIONS

- AMBIENT OPERATING TEMPERATURE:
 - SCAN ROOM: TEMPERATURE RANGE 64-79° F (18-26° C)
 - CONTROL ROOM: MAINTAIN TEMPERATURE AT 72° F (22° C)
 - EQUIPMENT ROOM (IF SEPARATE): TEMPERATURE RANGE 60-75° F (15-24° C)
- MAXIMUM TEMPERATURE RATE OF CHANGE OF 5° F (3° C)/HOUR.
- HUMIDITY: 30 TO 60 PERCENT NON-CONDENSING DURING OPERATION (ALL AREAS)
- MAXIMUM RELATIVE HUMIDITY RATE OF CHANGE IS 5 PER CENT RH/HOUR.
- ALTITUDE: NOT TO EXCEED 7875 FT. (2400M) ABOVE SEA LEVEL.
- THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
- DO NOT RESTRICT THE AIR INTAKE OR AIR EXHAUST OF THE SYSTEM COMPONENTS.
- ENVIRONMENTAL CONDITIONS LISTED ABOVE MUST BE MAINTAINED AT ALL TIMES INCLUDING FOR EXAMPLE OVERNIGHT, WEEKENDS, AND HOLIDAYS.

MAGNETIC INTERFERENCE SPECIFICATIONS

- CT Gantry must be located in ambient static magnetic fields of less than one Gauss to guarantee specified imaging performance. Ambient AC magnetic fields must be below 0.01 Gauss peak.
- CT computer equipment must be located in ambient static magnetic fields of less than ten Gauss to guarantee data integrity.
- CT control equipment must be located in ambient static magnetic fields of less than listed below to obtain specified geometric linearity.

CONSOLE/COMPUTER	10 GAUSS
CRT MONITOR	1 GAUSS
LCD MONITOR	50 GAUSS

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT LAYOUT
MODALITY TYPE: LIGHTSPEED RT16/XTRA
THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IT IS NOT TO BE USED FOR CONSTRUCTION PURPOSES. IT IS NOT TO BE USED FOR ACTUAL CONSTRUCTION PURPOSES; HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
6-64f
TYPICAL FINAL
(with GT1700 Table)


PROJECT: 6-64f
REVISION: 05
DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA

REVISION HISTORY:

SHEET
A1

TYPICAL WALL SUPPORT ELEVATIONS


S121



SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-0"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO STRUCTURAL DETAIL B7998D. NOTE: FOR SUPPORT LENGTHS LONGER THAN 96 IN. ADDITIONAL BRACING MAYBE REQUIRED.
2	THE PEDESTAL-CEILING MOUNT REQUIRES A FLUSH CEILING MOUNTING PLATE THAT IS STRUCTURALLY SUPPORTED TO HANDLE THE WEIGHT OF THE LOAD AS SHOWN IN DETAIL B50-31B. A MAVIG DESIGNED CEILING PLATE IS AVAILABLE AND CAN BE ORDERED - E9007WZ. THIS 14" x 14" PLATE HAS PRE DRILLED 1/2" HOLES IN EACH OF THE 4 MOUNTING CORNERS. IF AN EQUIVALENT PLATE IS USED, THE SUPPLIED TEMPLATE SHOULD BE USED TO DRILL THE REQUIRED 4 - 3/8" HOLES IN A PATTERN AS SHOWN IN DETAIL B50-31B. AN ADDITIONAL 3/8" HOLE IS REQUIRED FOR THE PEDESTAL-CEILING MOUNT SAFETY CHAIN.
3	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S121, FOR VARIAN POWER SUPPLY MODULE.
4	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE CEILING LASER POSITIONING LIGHT. SUPPORT TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL SQUARE AND IN THE SAME HORIZONTAL PLANE ABOVE FINISHED CEILING. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. SEE DETAIL B7998C.
5	MINIMUM 8'-0" x 2'-8" OPEN CEILING SPACE REQUIRED FOR OVERHEAD LASER INSTALLATION. CEILING TO BE FINISHED BY CUSTOMER/CONTRACTOR ONCE LASER INSTALLATION IS COMPLETE.
6	LEVELING AREA FOR GANTRY AND TABLE SEE DETAIL B78-17 ON SHEET S2.

CT - LIGHTSPEED RT16 XTRA Seismic Zone Anchoring Hardware
 (GANTRY) ANCHORS = Hilti KB3 - 5/8 x 10 in. (4 ea.)
 (GANTRY) Steel Inserts = P-11/16 x 1 x 2.5 in. (4 ea.)
 (TABLE) ANCHORS = Hilti KB3 - 5/8 x 10 in. (4 ea.)
 (TABLE) Steel Inserts = P-11/16 x 1 x 2.5 in. (4 ea.)
 (PDU) ANCHORS = Hilti KB3 - 3/8 x 3.75 in. (4 ea.)
 (CONTROL) ANCHORS = Hilti KB3 - 3/8 x 3.75 in. (4 ea.)
 ALL ANCHORS TO INCLUDE 1 FLATWASHER
 ALL BOLTS TO INCLUDE 6 FLATWASHERS, 1 LOCKWASHER AND 1 NUT
 ALL BRACKETS ARE SHIPPED WITH GE EQUIPMENT.

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED EQUIPMENT IS TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 1/4" BELOW THE FINISHED CEILING.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 1/8" IN 10'-0"
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.

SHEET TITLE: STRUCTURAL LAYOUT
 MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY ERRORS OR OMISSIONS IN THIS PLAN. THE USER OF THIS PLAN ASSUMES ALL LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
 6-64f
 TYPICAL FINAL
 (with GT1700 Table)

PROJECT	REVISION
6-64f	05

DATE: 25.Jan.12
 DRAWN BY: JGA
 CHECKED BY: JGA

REVISION HISTORY:


SHEET
 S1

GE Healthcare
 IS Services Design Center
 Milwaukee, Wisconsin

SUPPORT DETAIL
WALL SUPPORT FOR SIDE LASER

B7998D

REV. DATE: 03/15/03




DETAIL NOT TO SCALE

SUPPORT DETAIL
OVERHEAD COUNTERPOISED SUSPENSION

B50-31B

REV. DATE: 04/28/08




THE SUPPORT STRUCTURE FOR THIS CEILING MOUNTED OPTION AND A FLUSH MOUNTING PLATE MUST BE DESIGNED BY A STRUCTURAL ENGINEER AND INSTALLED BY A QUALIFIED CONTRACTOR PRIOR TO THE CT INSTALLATION

DRAWING NOT TO SCALE

SUPPORT DETAIL
TYPICAL CEILING SUPPORT FOR OVERHEAD LASER

B7998C

REV. DATE: 03/15/03




Note - Unless otherwise specified
1. This example is NOT recommended for concrete ceiling to false ceiling distance of more than 36" (914mm). Anti-sway devices will be required.

DETAIL NOT TO SCALE

CT GANTRY AND TABLE ANCHOR/LEVELING

B78-17

REV. DATE: 01/28/09



FLOOR VIBRATION

THE CT EQUIPMENT MAY BE SENSITIVE TO VIBRATION IN THE FREQUENCY RANGE OF 0.5 TO 20 HZ DEPENDING ON THE AMPLITUDE OF THE VIBRATION. IT IS THE CUSTOMERS RESPONSIBILITY TO CONTRACT A VIBRATION CONSULTANT OR QUALIFIED ENGINEER TO IMPLEMENT DESIGN MODIFICATIONS TO MEET THE SPECIFIC LIMITS. HOWEVER, IT IS ULTIMATELY THE CUSTOMER/ARCHITECT/ENGINEER RESPONSIBILITY TO DESIGN THE SITE SOLUTION.

STEADY STATE VIBRATION


THE MAXIMUM STEADY STATE VIBRATION TRANSMITTED THROUGH THE FLOOR SHOULD NOT EXCEED 10^{-3} M/S² RMS MAXIMUM SINGLE FREQUENCY ABOVE AMIENT BASELINE FROM 0.5 TO 80 HZ (MEASURED IN ANY 1 HOUR DURING A NORMAL OPERATING PERIOD).

TRANSIENT VIBRATION

THE BEHAVIORAL CHARACTERISTICS MUST BE SUCH THAT ANY MEASUREABLE TRANSIENT DISTURBANCE MUST ALSO BE MINIMIZED TO LESS THAN 0.01 M/S² PEAK-TO-PEAK.

EQUIPMENT LOCATION

TO MINIMIZE THE INTERFERENCE, THE SYSTEM SHOULD BE PLACED ON A SOLID FLOOR, LOCATED AS FAR AS POSSIBLE FROM THE VIBRATION SOURCES, SUCH AS PARKING LOTS, ROADWAYS, SUBWAYS, TRAINS, HALLWAYS, ELEVATORS, AND HOSPITAL PHYSICAL PLANTS. PLEASE NOTE THAT OTHER ITEMS NOT LISTED COULD ALSO BE POTENTIAL SOURCES OF VIBRATION.



DETAIL NOT TO SCALE

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. THE USER SHALL BE RESPONSIBLE FOR VERIFYING THE LOCATION OF ALL EQUIPMENT TO ACTUAL EQUIPMENT PURPOSES. HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

6-64f
TYPICAL FINAL
(with GT1700 Table)

PROJECT REVISION
6-64f 05

DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA

REVISION HISTORY:

SHEET


S2

SCALE: 1/4" = 1'-0"


ELECTRICAL PLAN

RECOMMENDED CEILING HEIGHT = 9'-0"


JUNCTION POINT DESCRIPTIONS




A




B




C



D



E



F

ELECTRICAL OUTLET LEGEND
CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS. HEIGHT ABOVE FLOOR DETERMINED BY LOCAL CODES UNLESS OTHERWISE SPECIFIED.

- ▲ DEDICATED TELEPHONE LINE(S) (SEE ELECTRICAL DETAIL ELEC-1 OR ELEC-67)
- ▲ NETWORK OUTLET (SEE ELECTRICAL DETAILS ELEC-83 AND ELEC-84 OR ELEC-87)
- ⊕ DUPLEX HOSPITAL GRADE, DEDICATED WALL OUTLET 120-V, SINGLE PHASE POWER
- ⊕ DUPLEX HOSPITAL GRADE, DEDICATED CEILING OUTLET 120-V, SINGLE PHASE POWER


DUCT HATCHING LEGEND

- ▨ ABOVE CEILING DUCT
- ▧ UNDER FLOOR DUCT
- ▩ TRENCH DUCT (FLUSH FLOOR)
- ▨ SURFACE FLOOR DUCT
- ▨ CABLE TRAY
- ABOVE CEILING CONDUIT
- BELOW FLOOR CONDUIT

JUNCTION POINT NOTES

- o ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER'S ELECTRICAL CONTRACTOR.
- o CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS
- o CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
- o CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
- o ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
 1. DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
 2. DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL POWER PURPOSES.
 3. DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
 4. PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
- o ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMER'S CONTRACTOR.
- o GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
- o 10 FOOT PITTAILS AT ALL JUNCTION POINTS.
- o ALL WIRING MUST BE THHN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT INSULATION. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- o GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.

PLEASE SEE BELOW FOR ADDITIONAL REQUIRED CONDUIT RUNS AND SIZES.



FEEDER TABLE - CT LightSpeed Pro 16 / RT / VCT

o CALCULATIONS BASED UPON NOMINAL VOLTAGE. WIRE SIZE IN AWG.
o RECOMMENDED FEEDER SIZES FROM DISTRIBUTION TRANS. TO POWER DISTRIBUTION UNIT.
o THE GROUNDING CONDUCTOR () WILL BE A 1/0 MINIMUM. THIS GROUND WILL RUN FROM THE EQUIPMENT BACK TO THE POWER SOURCE/MAIN GROUNDING POINT AND ALWAYS TRAVEL IN THE SAME CONDUIT WITH THE FEEDERS AND NEUTRAL.
o NEUTRAL MUST BE TERMINATED PRIOR TO OR INSIDE THE MAIN DISCONNECT PANEL AND NOT BROUGHT INTO THE POWER DISTRIBUTION UNIT.
o FOR A FULL SYSTEM UPS REFER TO ELECTRICAL DETAILS FOR UPS FEEDER WIRES.

RUN LENGTH IN FEET	POWER SUPPLY VOLTAGE					
	342-418 380	360-440 400	378-462 420	396-484 440	414-506 460	432-528 480
50	FEEDER GROUND	FEEDER GROUND	FEEDER GROUND	FEEDER GROUND	FEEDER GROUND	FEEDER GROUND
100	1/0 (1/0)	1/0 (1/0)	1/0 (1/0)	1 (1/0)	1 (1/0)	1 (1/0)
150	1/0 (1/0)	1/0 (1/0)	1/0 (1/0)	1 (1/0)	1 (1/0)	1 (1/0)
200	1/0 (1/0)	1/0 (1/0)	1/0 (1/0)	1 (1/0)	1 (1/0)	1 (1/0)
250	2/0 (1/0)	2/0 (1/0)	1/0 (1/0)	1/0 (1/0)	1 (1/0)	1 (1/0)
300	3/0 (1/0)	3/0 (1/0)	2/0 (1/0)	2/0 (1/0)	1/0 (1/0)	1/0 (1/0)
350	4/0 (1/0)	3/0 (1/0)	3/0 (1/0)	2/0 (1/0)	2/0 (1/0)	1/0 (1/0)
400	250M (1/0)	4/0 (1/0)	3/0 (1/0)	3/0 (1/0)	3/0 (1/0)	2/0 (1/0)

REV. DATE: 04/15/07

ADDITIONAL CONDUIT RUNS FOR ALL LIGHTSPEED AND BRIGHTSPEED SYSTEMS AND THE HISPEED QX/i (BY CONTRACTOR)

CONDUITS REQUIRED FOR INJECTOR: NEMOTO (CONDUITS ARE LOCATED ABOVE CEILING)

CONDUITS REQUIRED FOR BASE SYSTEM (CONDUITS ARE LOCATED ABOVE CEILING)

CONDUITS REQUIRED FOR UPS (CONDUITS ARE LOCATED ABOVE CEILING)

CONDUITS REQUIRED FOR LAP LASERLIGHTS (CONDUITS ARE LOCATED ABOVE CEILING)

CONDUITS REQUIRED FOR SMARTVIEW & SMARTSTEP MONITORS (CONDUITS ARE LOCATED ABOVE CEILING)

JUNCTION POINT DESCRIPTIONS

THE FOLLOWING MATERIALS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER'S ELECTRICAL CONTRACTOR

POINT	DESCRIPTION	QTY.	HARDWARE	DETAIL NO., SHT. E3
A1	MAIN DISCONNECT	1	PANEL - INCLUDED IN ORDER	ELEC-135
CL	PATIENT POSITIONING LASER LIGHT	1	COVERPLATE 4 X 4 X 4 IN. BOX 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
CLC	PATIENT POSITIONING LASER LIGHT CONTROL CONSOLE	1	COVERPLATE 4 X 4 X 4 IN. BOX 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8 ELEC-84
CTT	CT SCANNER	2	3 1/2 IN. DIA. BUSHING & LOCKNUT 2 1/2 IN. DIA. BUSHING & LOCKNUT	ELEC-9
ICC	INJECTOR CONTROL CONSOLE	1	EXTERNALLY CONNECTED	
IE	INJECTOR ELECTRONICS	1	6 X 6 X 4 IN. BOX 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD	1	SINGLE GANG BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
KP	KEYPAD	1	CONNECT EXTERNALLY TO 'CL'	
OC	OPERATORS CONSOLE	1	COVERPLATE 10 IN. DIA. CHASE NIPPLE 10 X 6 X 4 IN. BOX 2 1/2 IN. DIA. BUSHING & LOCKNUT 2 1/2 IN. DIA. BUSHING & LOCKNUT	ELEC-9
PM	POWER DISTRIBUTION UNIT	1	SPLIT COVERPLATE 3 1/2 IN. DIA. BUSHING & LOCKNUT 1 1/2 IN. DIA. BUSHING & LOCKNUT 16 IN. 90 DEGREE CONNECTOR 16 FT. LENGTH OF 1/2 IN. FLEXIBLE METAL CONDUIT SUITABLE CONNECTORS 12 X 16 X 4 IN. BOX 6 FT. LENGTH OF 1/2 IN. FLEXIBLE METAL CONDUIT	ELEC-9 ELEC-22
PSM	VARIAN POWER SUPPLY MODULE	1	VARIAN SUPPLIED, CUSTOMER/CONTRACTOR INSTALLED VARIAN POWER SUPPLY	
RGC	VARIAN RESPIRATORY GATING CAMERA	1	1 1/2 IN. DIA. CHASE NIPPLE	
RWS	VARIAN OPERATORS WORKSTATION	1	4 X 4 X 4 IN. BOX 1 IN. DIA. CHASE NIPPLE COVERPLATE	ELEC-8 ELEC-83 ELEC-84
SED	EMERGENCY OFF	1	SINGLE GANG 2 1/8 IN. DEEP FLUSH MOUNTED JUNCTION BOX.	ELEC-16
SVM	LCD MONITOR	1	COVERPLATE 6 X 6 X 4 IN. BOX 4 X 4 IN. OPENING IN COVERPLATE	ELEC-77
UPS	UPS CABINET	1	COVERPLATE 4 X 4 X 4 IN. BOX 1 IN. DIA. CHASE NIPPLE 1 1/2 IN. DIA. BUSHING & LOCKNUT (IF OPTIONAL 2 IN. CND IS USED, ADD THE FOLLOWING) COVERPLATE 4 X 4 X 4 IN. BOX 2 IN. DIA. BUSHING & LOCKNUT 2 IN. DIA. CHASE NIPPLE	ELEC-8
VED	VARIAN WALL SWITCH	1	SINGLE GANG 2 1/8 IN. DEEP FLUSH MOUNTED JUNCTION BOX.	ELEC-16
WL	WARNING LIGHT	1	X-RAY ON INCANDESCENT LIGHT FIXTURE DO NOT USE FLUORESCENT FIXTURES GE CAT. NO. WX1ABW-DF-XIU	
WLC	WARNING LIGHT CONTROLLER	1	E4502RL WARNING LIGHT CONTROL OR EQUIVALENT MAX 24V CONTROLLER	ELEC-72

CONTRACTOR SUPPLIED AND INSTALLED WIRING
ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS.

WIRE RUN, FROM - TO	QUANTITY, WIRE SIZE/COLOR
A1 > SED	1-ND 14 BLACK, 1-ND 14 WHITE, 1-ND 14 GREEN
A1 > PDU	3-BLACK, 1 GREEN - REFER TO FEEDER TABLE
480-V > A1	3 BLACK, 1 WHITE, 1 GREEN - REFER TO FEEDER TABLE
WLC > 1 PHASE	1-ND 14 BLACK, 1-ND 14 WHITE, 1-ND 14 GREEN
PM > WLC	1-ND 14 BLACK, 1-ND 14 WHITE, 1-ND 14 GREEN
WL > WLC	2-ND 14 BLACK, 1-ND 14 RED, 1-ND 14 WHITE

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL LAYOUT
MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. THE CONTRACTOR SHALL VERIFY ALL DIMENSIONS AND CONDITIONS OF THE SITE PRIOR TO ACTUAL EQUIPMENT INSTALLATION. GE HEALTHCARE SHALL NOT BE RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: 6-64f
TYPICAL FINAL
(with GT1700 Table)


PROJECT: 6-64f
REVISION: 05

DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA

REVISION HISTORY:

SHEET E1

INTERCONNECT DIAGRAM



POWER SPECIFICATIONS

CT LightSpeed Pro 16 / RT / VCT

(REV. DATE 04/15/07)
 VOLTAGE: PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS. RANGE OF LINE VOLTAGES: NOMINAL LINE VOLTAGE OF 380 TO 480, 3 PHASE, 50 OR 60 Hz.
 REQUIRED POWER SUPPLY: WYE-CONNECTED

MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

TABLE A ALLOWABLE INPUT VOLTAGES/CURRENT DEMAND

NOMINAL VOLTAGE	ABSOLUTE RANGE	CURRENT (AMPS)		MINIMUM STANDARD OVERCURRENT PROTECTION
		MAXIMUM	CONTINUOUS	
380	342-418	253	38	150-A
400	360-440	241	36	150-A
420	378-462	229	34	150-A
440	396-484	219	33	125-A
460	414-506	209	31	125-A
480	432-528	200	30	125-A

(ALL CALCULATIONS BASED UPON NOMINAL VOLTAGE)

PHASE-BALANCE: PHASE-TO-PHASE VOLTAGE MUST BE WITHIN +2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED LINE VOLTAGE AT A MAXIMUM DURATION OF 1 CYCLE AND FREQUENCY OF 10 TIMES PER HOUR.

VOLTAGE TRANSIENT OR IMPULSE ON THE INCOMING POWER MUST BE HELD TO A MINIMUM. TRANSIENTS CAUSED BY LIGHTNING, SURGES, LOAD SWITCHING, STATIC ELECTRICITY ETC. CAN CAUSE SCAN ABORTS OR, IN EXTREME INSTANCES, COMPONENT FAILURE IN THE COMPUTER SUBSYSTEM.

POWER DEMAND: CONTINUOUS POWER DEMAND = 25 KVA (MAX DEMAND = 150 KVA)

TABLE B MAXIMUM MOMENTARY POWER DEMAND.

DEMAND	CT HiSpeed
kVa *	150
POWER FACTOR AT	0.85

* DEMAND INCLUDES POWER FOR ENTIRE CT SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 6 PERCENT.

DISTRIBUTION TRANSFORMER: FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 225 KVA. GE DOES NOT RECOMMEND USING A REGULATION DEVICE.

NOTE: THE CT SYSTEM MUST NOT BE POWERED IN A MULTIPLE INSTALLATION WHERE FILM CHANGERS ARE USED. FILM CHANGERS UTILIZE A LARGE NUMBER OF HIGH POWERED CLOSELY SPACED EXPOSURES WHICH MAY COINCIDE WITH THE CT SCAN.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

DIAGRAM KEY

---	CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN ADEQUATE CONDUIT OR RACEWAY.
---	GE FURNISHED CABLE RUNS. ROUTE IN EMPTY CONDUIT OR RACEWAY.
59' [18M]	MAXIMUM RUN LENGTH BETWEEN JUNCTION POINTS. Feet [Meters]

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

GE Healthcare
 IS Services Design Center
 Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL SPECIFICATIONS
 MODALITY TYPE: LIGHTSPEED RT16/XTRA
 THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. THE USER SHALL VERIFY THE LOCATION OF ALL EQUIPMENT AND THE LOCATION OF ALL ELECTRICAL EQUIPMENT BEFORE INSTALLATION. THE USER SHALL BE RESPONSIBLE FOR OBTAINING ALL NECESSARY PERMITS AND FOR OBTAINING ALL NECESSARY CONSTRUCTION PERMITS. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
 6-64f
 TYPICAL FINAL
 (with GT1700 Table)

PROJECT	REVISION
6-64f	05

DATE: 25.Jan.12
 DRAWN BY: JGA
 CHECKED BY: JGA


REVISION HISTORY:

SHEET
 E2

ELECTRICAL DETAIL
BOX WITH COVERPLATE (TYPICAL)

ELEC-8

REV. DATE: 09/30/94




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
EMERGENCY OFF BUTTON

ELEC-16

REV. DATE: 08/22/05




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH COVERPLATE AND NETWORK JACK

ELEC-83

REV. DATE: 10/06/98




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
NETWORK CONNECTION (TYPICAL)

ELEC-84

REV. DATE: 03/06/04




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
WARNING LIGHT DIAGRAM

ELEC-72

REV. DATE: 08/12/08




UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER, ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

DRAWING NOT TO SCALE

ELECTRICAL DETAIL
CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9

REV. DATE: 08/08/94




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH SPLIT COVERPLATE (TYPICAL)

ELEC-22

REV. DATE: 10/13/94




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
INSITE CONNECTION (TYPICAL)

ELEC-1

REV. DATE: 04/24/02




ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
MAIN DISCONNECT PANEL

ELEC-135

REV. DATE: 09/27/04




TYPICAL CONDUIT AREA

EQUIPMENT DETAIL
CEILING PLATE (TYPICAL)

ELEC-77

REV. DATE: 08/13/07



DETAIL NOT TO SCALE

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. THE CONTRACTOR SHALL BE RESPONSIBLE FOR VERIFYING THE LOCATION OF ALL UTILITIES TO ACTUAL EQUIPMENT SCHEDULED TO BE INSTALLED. GE IS NOT TO BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
6-64f
TYPICAL FINAL
(with GT1700 Table)

PROJECT	REVISION
6-64f	05

DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA


REVISION HISTORY:

SHEET
E3

EQUIPMENT DETAIL
SIDE MOVING LASER

B7998B

REV. DATE: 03/15/03




DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INJECTOR ON OVERHEAD
COUNTERPOISED SUSPENSION

B50-31

REV. DATE: 05/04/02

VERIFY MOUNTING ASSEMBLY DIMENSIONS
WITH INJECTOR MANUFACTURER.




DRAWING NOT TO SCALE

EQUIPMENT DETAIL
NEMOTO INJECTOR ELECTRONICS

B8007C

REV. 00: 7/9/08




DISPLAY UNIT DRAWING NOT TO SCALE MAIN UNIT

EQUIPMENT DETAIL
VARIAN POWER SUPPLY (PSM)

R4504D

REV. DATE: 07/14/08




DETAIL NOT TO SCALE

EQUIPMENT DETAIL
GATING INSTALLATION - TYPICAL ROOM SECTION

E8819KA

REV. DATE: 01/23/09




SIDE VIEW
DETAIL NOT TO SCALE

EQUIPMENT DETAIL
OVERHEAD LASER

B7998A

REV. DATE: 03/15/03




DETAIL NOT TO SCALE

EQUIPMENT DETAIL
LAP LASER KEYPAD

B7998F

REV. DATE: 05/23/03




DETAIL NOT TO SCALE

EQUIPMENT DETAIL
10KVA UPS SYSTEM

B79-99ZA

REV. DATE: 04/18/06




UPS CABINET

EQUIPMENT DETAIL
POWER DISTRIBUTION UNIT

B78-58D

REV. DATE: 01/28/09




FRONT VIEW

SIDE VIEW

EQUIPMENT DETAIL
OPERATORS CONSOLE

B78-58A

REV. DATE: 01/28/09




FRONT VIEW

SIDE VIEW

EQUIPMENT DETAIL
MAIN LINE CONTACTOR - OPTIONAL

E45-02AE

REV. DATE: 02/20/04



FRONT VIEW


SIDE VIEW

BACK VIEW

EQUIPMENT DETAIL
MONITOR ON OVERHEAD
COUNTERPOISED SUSPENSION

B5031S

REV. DATE: 04/28/08



DRAWING NOT TO SCALE

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. THE USER SHALL BE RESPONSIBLE FOR VERIFYING THE LOCATION OF THE EQUIPMENT TO ACTUAL CONSTRUCTION PURPOSES; HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.










PROJECT TITLE:
6-64f
TYPICAL FINAL
(with GT1700 Table)

PROJECT REVISION
6-64f 05
DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA

REVISION HISTORY:

SHEET
D1

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin



GE Healthcare

IS Services Design Center

Wisconsin
Milwaukee

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: LIGHTSPEED RT16/XTRA

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. GE HEALTHCARE SHALL NOT BE RESPONSIBLE FOR THE DESIGN OF THE BUILDING OR FOR THE ACTUAL CONSTRUCTION OF THE EQUIPMENT. GE HEALTHCARE SHALL NOT BE RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

6-64f

TYPICAL FINAL

(with GT1700 Table)

PROJECT	REVISION
6-64f	05

DATE: 25.Jan.12
DRAWN BY: JGA
CHECKED BY: JGA

REVISION HISTORY:

SHEET

D2