

Definium™ XR120

See digitally. See beyond.



Definium™ XR120 brings the world of digital technology to imaging facilities that rely on analog equipment. As a competitive digital alternative, the Definium XR120 wireless digital detector seamlessly adapts to both fixed and mobile x-ray systems and is available with the latest 127 μ m cesium iodide technologies.

Advantages:

- Anywhere, anytime access
- Versatility starts with flexibility
Available for both fixed and mobile devices
- Fast preview
Image acquisition within two seconds
- Intuitive interface
- Available in multiple sizes
10" x 12" | 14" x 17" | 17" x 17"
- Advanced stitching capabilities
- Chiropractic tools are built right in
endless imaging capabilities with Smart Guide tutorials
- Durability you can depend on
Water-resistant to most typical water spills and tested for dust resistance (Ip66)
- Detector payload and drop grade
support up to 300 kg (661 lb) over the whole area of detector window and 150 kg (330 lb) partial pressure. It can also withstand drops from heights of up to 1.2 meters
- Printing made easy
- More mobility
A lightweight design* with wireless detectors
- User-friendly design
curved edges and an anti-slip grip nonslip back



Request a call back